

Ultraviolet Radiation Device Use for Decontamination
Testing Trial Report

Tested 3 discharged patient rooms with the UV device after the room was cleaned by EVS.
 Tested the rooms with the Hygiena ATP swab.
 Set up the UV device I in the patient room and 1 in bathroom with the appropriate time for bacteria.
 Ater the dose of UV radiation is completed, the room was tested with the Hygiena ATP swab.
 The results showed an improvement on surfaces tested

High Touch Surfaces	ATP Test Results	After UV Device ATP Test Results
Bed Rail Controls	9	3
	12	5
	1	
Tray Table	17	2
	26	12
	20	16
IV Pole (grab area)	12	12
	2	1
	7	10
Telephone	1	2
	8	5
	5	1